

# ENTEROCOCCI (10 TESTS)

Product number: 700285



Enterococci tests are a proven indicator for assessing microbiological water quality. The presence of enterococci in potable or ballast water may signal fecal contamination and potential hygienic risk—making regular testing critical for compliance, crew health, and operational confidence. To ensure accurate and consistent results, samples must be incubated at the correct temperature using a reliable incubator.

**Test Parameters** - Developed to detect *Enterococcus* species in potable and ballast water samples, delivering clear results within **24 hours**.

**Quantity** - Each kit contains **10 tests per box**, providing excellent value for routine onboard testing.

**What's Included:** the kit is supplied with 10 x 100ml Thio bags. Everything required to run multiple tests is conveniently packed in one complete solution.

**Protect your water systems. Support compliance. Safeguard crew wellbeing.**

Choose a simple, effective testing solution you can rely on—on every voyage.

## Instructions for use

### Step 1 – Sampling and incubation

Add reagent to sample and incubate at 35°C ± 0.5°C for 18 hours with WSS Portable Incubator (777930)

### Step 2 - Results

Read results:

Blue = Negative

Green = Positive – Enterococci Present

## To be used with

| Product                           | Product Number |
|-----------------------------------|----------------|
| Potable Water Test Kit W/O E.COLI | 777619         |
| Ballast Check 2 Test Kit          | 777848         |

**Wilhelmsen Ships Service**  
Phone: (+47) 67 58 40 00

Postal Address:  
PO Box 33, NO-1324  
Lysaker,  
Norway  
wilhelmsen.com

V1.0-042026

Contact your local WSS customer service for prices and availability

Scan now



for a closer look